

User guidelines for the T2*- estimation MATLAB Tool

General information

The tool consists of **2 components**:

1. temporal voxel-by-voxel noise filtering of 3 gradient echoes acquired as 3-D multi-volume fMRI data
2. voxel-by-voxel estimation of T2* at all acquired repetition times

These 2 steps are **not independent**: Step 2, the T2* estimation, can not be done before the input data have been processed by Step 1, the temporal voxel-by-voxel noise filtering.

The **input data** must contain three .nii files named *Echo1*.nii, *Echo2*.nii and *Echo3*.nii files containing fMRI echo series, as well as three .mat files named *Echo1*.mat, *Echo2*.mat and *Echo3*.mat containing DICOM information related to *Echo1*.nii, *Echo2*.nii and *Echo3*.nii.

Location of the datasets: the datasets used so far are placed on the CBIA **gryf** data server in the directory
/data/common/PROJECTS/_Current/CEITEC_T2star/T2star_data.

An **example dataset** that can be immediately used by the CBIA as input to temporal voxel-by-voxel noise filtering is located at
/data/common/PROJECTS/_Current/CEITEC_T2star/T2star_data/Example

Installation of the tool: simply place the T2star_github directory to a location of your choice.

Startup of the tool:

- launch MATLAB >=2016b
- change to your T2star_github directory

The tool was developed and tested using **MATLAB 2016b** under **Windows 7** and **Windows 10**.

1. Voxel-by-voxel noise filtering of the time-course of 3 fMRI gradient echoes.

This part of the tool is described in a paper to appear in ...

1.1 the noise problem

...

1.2 *The denoising algorithm.*







...

1.3 *Running the MATLAB implementation of the denoising algorithm.*

...

The denoising algorithm is encapsulated in the script named T2star_TV_T2star github.m which calls several other scripts and one mex-function.

On completion of each echo's denoising, the algorithm creates a summary plot which helps pinpoint parameters that should be adjusted to obtain optimum convergence, as explained further on:



The denoising script is launched by the command

```
>> T2star_TV_github
```

after which a self-explanatory dialog follows.

Algorithm notifications are usually introduced by "CBIA>>" to distinguish them from system messages.

CBIA>>please make sure spm12 toolbox is on your MATLAB path

Prompts for user actions are mostly terminated by a colon:

CBIA>>please enter full directory path containing

xxx_Echo_N.nii and

yyy_Echo_N_dicom_header.mat,

such as I:\Data\T2star\4055B\FCNI1_mag\cwu:

while informative messages about the algorithm progress begin with an ellipsis (...):

CBIA>>...unfiltered echoes will be read from

I:\Data\T2star\4055B\FCNI1_mag\cwu

CBIA>>...reading 4-D echoes Echo_1, Echo_2, Echo_3

and Echo_1_dicom_header, Echo_2_dicom_header, Echo_3_dicom_header

CBIA>>...reading TE(1),TE(2),TE(3) and TR from dicom headers

CBIA>>...concatenating 4-D echoes Echo_1, Echo_2, Echo_3

to 5-D datacube Echoes with echo number as the 5-th dimension

Denoising is performed on 3-D volumes of all 3 echo trains for all fMRI-acquired time instants.

A z-slice is selected here for 2-D visual representation only, and does not affect denoising:

CBIA>>pick a z-slice, or press RETURN for default: [23]:


To speed up computations by parallel processing on multi-core CPUs, the number of threads of the user's CPU can be entered:

CBIA>>enter number of CPU threads, or press RETURN for 8 threads:

Setting the algorithm parameters:


CBIA>>enter $\log_2(\mu)$, or press RETURN for default [-16]:

The parameter $\mu=2^{\log_2(\mu)}$ defines weighting of the L2-norm (mean square error) between the measured and the denoised echo time-course: the smaller μ , the smoother, in terms of total variation **TV**, will the denoised signal be. μ affects the speed and the stability of the TV optimization process.



CBIA>>enter `log2(beta)`, or press RETURN for default [-6]:

The parameter `beta=2^log2(beta)` affects the speed with which the true total variation TV will approach its approximation `w_norm`: the larger `beta`, the faster will TV and `w_norm` meet. `beta` affects the speed and the stability of the TV optimization process.




CBIA>>enter nr_shrinks, or press RETURN for default [10]:

The parameter `nr_shrinks` defines the number of **shrinkage** iterations, i.e. those in which total variation TV is reduced. `nr_shrinks` affects the accuracy of the TV optimization process.


CBIA>>enter nr_grads, or press RETURN for default [10]:

The parameter `nr_grads` defines the number of **steepest descent** operations, i.e. those in which square error `data_term` between the filtered and the noisy echo is reduced. `nr_grads` affects the accuracy and the stability of the `data_term` optimization process.

The product `nr_shrinks*nr_grads` should be large enough for the `data_term` to reach the plateau:



Also, the product `nr_shrinks*nr_grads` should be large enough for the **primal** and the **dual** functional to meet:



```
-----  
CBIA>>...TV filtering results will be stored in  
I:\Data\T2star\4055B\FCNI1_mag\cwu\TV_filtered_echoes\ -16  -6  10  10  
-----
```

The algorithm progress is reported:

```
CBIA>>...temporal TV filtering of Echo_1  
shrink_count=1  
shrink_count=2  
shrink_count=3  
shrink_count=4  
shrink_count=5  
shrink_count=6  
shrink_count=7  
shrink_count=8  
shrink_count=9  
shrink_count=10  
ElapsedTime= 176.0191  
CBIA>>...temporal TV filtering of Echo_2  
shrink_count=1  
shrink_count=2  
shrink_count=3  
shrink_count=4  
shrink_count=5  
shrink_count=6  
shrink_count=7  
shrink_count=8  
shrink_count=9  
shrink_count=10  
ElapsedTime= 176.7533  
CBIA>>...temporal TV filtering of Echo_3
```




```
shrink_count=1
shrink_count=2
shrink_count=3
shrink_count=4
shrink_count=5
shrink_count=6
shrink_count=7
shrink_count=8
shrink_count=9
shrink_count=10
ElapsedTime= 178.1268
CBIA>>...saving TV filtering results to
      I:\Data\T2star\4055B\FCNI1_mag\cwu\TV_filtered_echoes\ -16  -6  10  10
-----
```

after which the termination notification appears:

```
CBIA>>...temporal TV filtering of the 5D datacube has been completed.
```

Finally the user is prompted to peruse the denoising results:

To visually compare unfiltered and TV-filtered echoes at voxels of your choice,
pick repeatedly those you are interested in.



```
pick_another_pixel:[return/0]
pick_another_pixel:[return/0] 0
>>
>>
```

2. voxel-by-voxel estimation of T2* at all acquired repetition times.

2.1 Estimation of T2* by weighted exponential matching of the three TV-filtered gradient echoes.

The MATLAB tool is based on the original implementation of the algorithm described in *Michálek et al.*: “Fast and accurate compensation of signal offset for T2 mapping”, <https://link.springer.com/article/10.1007/s10334-019-00737-3>, which processed only single 2D slices of the input MRI data. For use with 3D fMRI data, the original code has been enhanced to include also the z-dimension (volume) and *t*-dimension (time). Resulting T2* estimates and noise masks are output as 4-D .nii files.

2.2 Running the T2*estimation MATLAB script.

The T2*estimation script is launched by the command

```
>> T2star_series_WLS_MFL
```

after which a self-explanatory dialog consisting of informative messages (introduced with ...) or prompts for user action (terminated with:) follows:

```
CBIA>>...please make sure spm12 toolbox is on your MATLAB path
CBIA>>please enter full directory path containing TV-filtered echoes
      TV_Echo_1.nii, TV_Echo_2.nii and TV_Echo_3.nii
      such as I:\Data\T2star\4055B\FCNI1_mag\cwu\TV_filtered_echoes\ -16   -6   10   10:
CBIA>>...TV-filtered echoes will be read from
      I:\Data\T2star\4055B\FCNI1_mag\cwu\TV_filtered_echoes\ -16   -6   10   10
CBIA>>...reading 4-D TV-filtered echoes TV_Echo_1, TV_Echo_2, TV_Echo_3
CBIA>>...concatenating 4-D TV-filtered echoes TV_Echo_1, TV_Echo_2, TV_Echo_3
      to 5-D datacube Echoes with echo number as the 5-th dimension
CBIA>>...unfiltered echoes will be read from
      I:\Data\T2star\4055B\FCNI1_mag\cwu
CBIA>>...reading 4-D unfiltered echoes Echo_1, Echo_2, Echo_3
      and Echo_1_dicom_header, Echo_2_dicom_header, Echo_3_dicom_header
CBIA>>...concatenating 4-D unfiltered echoes Echo_1, Echo_2, Echo_3
      to 5-D datacube Ref_echoes with echo number as the 5-th dimension
CBIA>>...reading TE(1),TE(2),TE(3) and TR from dicom headers
CBIA>>maximum slice:46
```


T2* estimation is performed on whole 3-D volumes for all fMRI acquired time instants.

z-slice is selected here for 2-D visual representation only, and does not affect denoising:

```
CBIA>>pick a z-slice for 2D visualizations, or press RETURN for default: [23]:  
slc=23
```

Before starting the T2* estimation, a gray-value **threshold** separating background pixels from valid fMRI information is generated automatically. The threshold generation relies on the assumption that the histogram of Echo_1 array (i.e. the brightest echo) exhibits two distinct peaks. The user is presented a plot with the histograms as well as the threshold-generated background mask, and is given the possibility to change the threshold:

```
CBIA>>press RETURN to accept the threshold, or enter a different one: [12692.5]:
```



```
-----  
CBIA>>...interleaving 4-D TV-filtered echoes TV_Echo_1, TV_Echo_2, TV_Echo_3  
        according to their absolute acquisition time  
        to 4-D datacube interleaved_echoes  
-----
```

```
CBIA>>...calculating T2star values for all volumes of input data  
-----
```

The total execution time needed to estimate T2* for all voxels at all measured fMRI volumes is reported
ElapsedTime=34.8382

and the directory where the 4-D results are saved as nifti files is displayed:



```
CBIA>>...saving noise mask and unmasked/masked T2star values as 4-D .nii files to  
        I:\Data\T2star\4055B\FCNI1_mag\cwu\TV_filtered_echoes\ -16  -6  10  10\T2star_series_WLS_MFL  
-----
```

For the user-selected slice, some raw/limited/masked 3-D T2* maps and a noise mask are generated, visualized and saved as multipage tiffs in a slice-related directory. These can later be viewed, e.g., using *ImageJ*:

```
CBIA>>...extracting 3-D data at your selected slice number 23 for visualization  
masked_pixel_ratio = 0.41266  
CBIA>>...saving raw/limited/masked 3-D T2stars and noise mask at your selected slice number in  
I:\Data\T2star\4055B\FCNI1_mag\cwu\TV_filtered_echoes\ -16  -6  10  10\T2star_series_WLS_MFL\slc23  
-----
```

And finally the user is given the possibility to inspect the T2* time courses at voxels of their choice:

```
CBIA>>...To visually inspect echo and T2star time courses at voxels of your choice,  
        pick repeatedly those you are interested in.
```



```
pick_another_pixel:[press RETURN for YES/0 for NO]
pick_another_pixel:[press RETURN for YES/0 for NO]
pick_another_pixel:[press RETURN for YES/0 for NO] 0
```

>>